




CME SOIL SAMPLING WITH WIRELINE

The **CME wireline punch core system** allows you to collect soil samples ahead of the core bit while using a core barrel. For sampling sandy or silty soils, sticky clays or compacted shales where contact by the drilling fluid is acceptable contamination, a plastic basket retainer is used to retain the sample. Move from soil sampling to rock coring without tripping wireline drill rods.



Available for all N, H & P size wireline and conventional core barrels. Tungsten carbide bit and 1.5" punch shoe shown. Other core bit types and punch shoe lengths available.

[CLICK HERE FOR PART #'s](#)



For more information, call 800-325-8827
or email productupdates@cmeco.com!

DEDICATION QUALITY KNOWLEDGE

CME is dedicated to providing the highest quality products and superior customer service, supported by the most experienced and knowledgeable sales staff in the industry.